

## FIXATION



**PRE-OP**



**FIXATION DRILL GUIDE**




**SEAT FIXATION GUIDE**

## S.I.M. SEATING



**CUSTOM S.I.M. DEVICE**



**S.I.M. DEVICE SEATED**




**IO SCAN WITH S.I.M. DEVICE**

## SURGERY



**GUIDED OR  
FREE-HAND SURGERY**

**SEAT MUAS AND  
SUTURE PATIENT**



**PLACE TI-BASE FREE SCAN FLAGS**  
CALL A NAVIGATION REPRESENTATIVE  
FOR A COMPLETE LIST OF OPTIONS!



**IO SCAN S.I.M. DEVICE  
AND SCAN FLAGS**

## DELIVERY



**STL DESIGN MATCH**



**PRINTED STL DESIGN**



**TRANSITIONAL FP1 DELIVERY**



# EZ Load™ CONVERSION

SIMPLY PRINT & DELIVER! NO ACRYLIC, ADJUSTING, OR BITE CHECKS



## WHAT IS EZ LOAD™?

EZ Load™ is a proprietary, standalone conversion process that eliminates the need for the traditional, arduous chairside processing of the conversion prosthesis after surgery. No luting or bite checks, deliver the prosthesis, and release the patient. EZ Load™ can be used in combination with free-hand or fully guided surgical procedures. This system enables our team to pre-design the conversion prosthesis in a fully digital workflow, and data archiving creates a disposable prosthesis – no more repair of a broken conversion, simply re-print and replace.

After placement of the MUA's, the surgeon sutures and scans the appropriate scan flags with an intra-oral scanner (no special camera or technical chairside assistance required). The proprietary S.I.M. (Smart Implant Mapping) device creates an accurate and fast IO data capture. The NavaGation team, through model matching, will adapt the diagnostic design to the implant positions and either email the STL file for in-office printing or fabricate the prosthesis in the lab for later delivery.

The design process is based on the patient's pre-op VDO (or adjusted bite as directed by the clinical team), and all of this is accomplished during the planning phase, simplifying the conversion process and greatly shortening the surgery. This precise, digital, diagnostically driven, conversion prosthesis serves as a healing device and doubles as a first-stage prototype that drives the final smile design.

Data archiving enables the restorative team to achieve a flawless final zirconia hybrid restoration in as little as two clinical appointments.

## WHAT WE NEED

### EDENTULOUS

- ✓ PVS Reline Denture
- ✓ 360° STL Denture Scan
- ✓ Opposing Arch and Bite Scans
- ✓ CBCT

### DENTATE

- ✓ Full-Arch Pre-Op IO Scans
- ✓ CBCT

### FULL-FACE SMILE PICTURES

- ✓ Resting
- ✓ Smiling
- ✓ Retracted
- ✓ Profile

**FREEHAND SURGERY DATA MUST BE SUBMITTED 10 DAYS PRIOR**  
**GUIDED SURGERY DATA MUST BE SUBMITTED 4 WEEKS PRIOR**

## DAY OF SURGERY

- ✓ Submit data through portal at **navagation.net**
- ✓ **Do not** direct upload through scanner portal (e.g., TRIOS, Primescan, etc.)
- ✓ Call us at **(844)-628-2428** or text your NavaGation rep. to confirm data upload

## DELIVERY OPTIONS

- ✓ STL Transfer In-Office Printing  
FP1 - 30 to 35 minute design turnaround  
FP2 or FP3 - 45 to 50 minute design turnaround  
*\* estimated times per arch*
- ✓ Next-Day Delivery Within 24 Hours